


# Introducing LiDCOview<sup>SE</sup>



NEW

The easy way to review LiDCOplus and LiDCOrapid data:

LiDCOview Standard Edition (SE) is designed as a graphical display of historical LiDCOplus and LiDCOrapid hemodynamic data. LiDCOview<sup>SE</sup> provides the clinical researcher with easy to view beat-to-beat hemodynamic data that was collected with the LiDCOplus and LiDCOrapid.



## Key Benefits:

- Research
- Audit
- Training
- Case Study

Download  
LiDCOview<sup>SE</sup>  
Free  
[www.lidco.com](http://www.lidco.com)

## Key Features of LiDCOview<sup>SE</sup>

- Easy to view graphic display of hemodynamic data
- Simple Navigation: zoom from 2min to 24hr with Calibration, Event and CVP flags displayed
- Data selection for numeric values and analysis
- Numeric values - includes dynamic preload response values: SVV and PPV
- Printing or screen shot for inclusion in papers and presentations
- On screen help
- Indexed values and real-time stamp with LiDCOplus Version 4 software



LiDCOview<sup>SE</sup> is intended for review of historical data for research and education.  
It is **not** intended for use in the diagnosis of disease or other conditions.

## Navigation Controls:

Overview and Zoom Bars make finding data easy.

**Yellow Highlighted Section indicates which portion of data is displayed on the Zoom Control bar; Up to 24hours can be selected**

**Overview Control Bar: 24hours of data with event flags displayed**

**Time is displayed both in absolute and relative terms**

**Flags are replicated for ease of navigation**

**Zoom Control Bar (Yellow): Select either end to adjust zoom; select middle to move**

**Press to advance by 2hrs or 24hrs**

Overview: Up to 24hrs in 2hr increments - Total file length: 23:53:27

Zoom: Enlarge, reduce or move - Range Selected: 02:19:23

Legend: No Data, CVP/RAP Value, Event Number, Calibration Factor

Hemodynamic Data is displayed as a Trend.  
Highlight a selection to see the specific numeric values.

**Flags correspond between Zoom Bar and Trend graphs**

**Highlighted section for numeric values; duration of selection**

**Y-axis adjustment**

**X-axis (time)**

**Numeric values**

Zoom: Enlarge, reduce or move - Range Selected: 02:19:23

Highlighted Section: 00:03:37

DO<sub>2</sub> ml min<sup>-1</sup> 588

CO l min<sup>-1</sup> 4.0

LiDCOview<sup>SE</sup> will run on the following Windows operating systems: NT, 2000sp3, XPsp2 and in 32 bit mode in Vista and 7. It requires Windows Installer and Internet Explorer 5.1. LiDCOview<sup>SE</sup> runs inside the Windows .Net 2.0 run-time environment, which requires 280Mb of hard disk space. LiDCOview<sup>SE</sup> uses the heart beat data file (\*.HBD, Version 3) or the lidcoview (\*.LVU, Version 4) available from the LiDCOplus or LiDCOrapid Engineering Screen.

## LiDCOview<sup>PRO</sup>

with even more features including:

- Advanced Navigation and Information Display
- Display: LiDCO Curves, Cal and Event History
- Analysis and Reporting
- Multiple File Display
- Graphics and Excel Exporting

Also Available

Email via: [www.lidco.com](http://www.lidco.com)  
Global Contact: +44 (0) 1223 830666